

First Position




Seating Concepts™

NEW

Balanced Aire Wheelchair Cushion™

Adjustable Skin Protection Cushion

The Balanced Aire Wheelchair Cushion™ is our Premier Pressure management wheelchair cushion. The Balanced Aire Wheelchair Cushion™ offers you dynamic and aggressive positioning, unmatched stability and a highly advanced level of prevention and treatment of pressure ulcers.



Mason Medical Products


Division of MRC Industries, Inc.

Tel: 800/233-4454 • 516/328-6900 • Fax: 516/328-6622

www.masonmedical.com

“Balanced Aire Wheelchair Cushion™”

Adjustable Skin Protection Cushion



FEATURES AND BENEFITS

- Flexible interconnected, independent 4" air cells allow for deep immersion into the cushion providing unvarying low pressure and an increased blood supply necessary for healthy tissue maintenance.
- The cells shift and shape to the seat area providing a low friction, low shear surface.
- The waterproof stretch cover is breathable providing a cool comfortable seating environment.
- The cushion is made of neoprene which is flame resistant and latex free.
- The easy to use pump allows for exact fill rate for maximum comfort and support.
- Weight capacity 350lbs
- 2 year non-prorated warranty

ITEM #:	8047-16	16" x 16" x 4"
	8047-18	18" x 16" x 4"
	8047-20	20" x 16" x 4"

Low Profile Cushions

ITEM #:	8047-16-2	16" x 16" x 2"
	8047-18-2	18" x 16" x 2"
	8047-20-2	20" x 16" x 2"



Mason Medical Products

Division of MRC Industries, Inc.

Tel: 800/233-4454 • 516/328-6900 • Fax: 516/328-6622