



PHYSICAL THERAPISTS RECOMMEND HOME SAFETY PRODUCTS

Accredited Medical Equipment & Supply



Product # INV9670

I-Class™ Blow-Molded Transfer Bench

- Make transfers in and out of the bathtub safe and comfortable.
- Extra sturdy, blow-molded bench holds up to 400 lbs.
- Reversible backrest snaps in and out to accommodate right or left transfers.
- Textured, nonslip seat has drain holes and built-in slots for hand-held showers and soap.
- Durable 1" aluminum frame will not rust.

Quad Canes

Small Base #ISG10110A

Large Base #ISG10111A

Handle assemblies may be rotated for left- or right-hand use. Stable chrome-plated steel base; offset handle with orthopedic handgrip provides better fit and added comfort. Height adjusts in 1" increments from 31" - 36"H. Weight capacity: 250 lbs.



Adult Hemi Walker
PRODUCT #CEXA873CO

This lightweight walker combines support and stability. Easy to lift and maneuver, even with one hand. Extra handgrip provides ease in folding. Two-level, comfortable handgrips provide support for standing from a seated position. Slip resistant rubber tips. Height adjusts from 29(1/2)" - 34(1/2)". Depth when folded is 3". Weight: 2.85 lbs. Weight capacity: 250 lbs.

OTHER ITEMS TO CONSIDER

- LONG HANDLE BRUSHES
- LONG-HANDLED SHOE HORN
- ELASTIC SHOE LACES
- REACHERS
- SLIDING BOARDS
- MALE & FEMALE URINALS
- HAND HELD SHOWER HEAD
- COMMUNE
- STRAIGHT CANE
- WALKER
- GAIT BELT

(800) 974 1234 Ext. 551

www.accreditedmedicalequipment.com

CONSIDER THESE ITEMS . . .



Medline Rollator The Walker You Can Ride On

PRODUCT #MDS86810BLK

This rolling "walker" has a padded seat, curved backrest, and loop brakes for quick and easy stopping. It's easy to maneuver and folds for storage. Hand grips adjust to your height: measure 32" from floor in lowest position or 36" from floor in highest position. Pull on grips to stop, push to lock brakes. Seat is 14" x 14"; 22"H from floor. When open, unit measures: 24"W x 22"L and weighs 16 lb. Capacity: 250 lb. Basket comes standard.



Handheld Shower PRODUCT #MDS41000W

Patient Therapy/Rehabilitation > Aids to Daily Living > Bathing Aids

Handheld Shower Allows Users to Sit During Showers

Equipped with a 60" hose (an extension hose).

Allows users to take a shower while sitting on a shower seat or transfer bench.

Simple installation.



Chrome Grab Bars PRODUCT #MDS86012CHR

Patient Therapy/Rehabilitation > Bath Safety & Commodes > Bathroom Grab Bars

Chrome Grab Bars Help Prevent Falls in the Shower

Essential for the elderly or wheelchair-bound, shower bars are a great safety measure for any shower. 1 1/2" knurled chrome diameter, chrome-finished bars may be used horizontally or vertically. Complete with mounting hardware which includes two rubber pads and six 2" rust resistant screws per bar.

(800) 974-1234 Ext. 551

www.accreditedmedicalequipment.com