


Multi-Ply™ Series 6500

VARIABLE DENSITY SLEEP SYSTEM

Multi-Ply 6500

6500-1 76"
6500-2 80"
6500-3 84"


6yr
warranty



Multi-Ply 6500-Lite

6500-LT 1 76"
6500-LT 2 80"
6500-LT 3 84"


4yr
warranty



Multi-Ply 6500-Global

6500-GL 1 76"
6500-GL 2 80"
6500-GL 3 84"

15yr
warranty



* Our Multi-Ply 6500 mattress line is of the same construction as mattresses Mason Medical previously designed and manufactured that were sold exclusively by Direct Supply, Inc., under Direct Supply's Panacea® mark prior to December 31, 2004, The Multi-Ply 6500 used to be sold as The Panacea®; the Multi-Ply 6500-Lite used to be sold as The Panacea® Tender; and the Multi-Ply 6500-Global used to be sold as The Panacea® Comfort Max. Mason Medical no longer manufactures Panacea® brand mattresses for Direct Supply. Our Multi-Ply 6500 mattresses are not approved, endorsed, or licensed by Direct Supply.

Mason Medical's exclusive manufacturing agreement with Direct Supply for the Panacea® Mattress line has expired. These Mason Medical manufactured mattresses are now available through all other distributors.

All Mattresses are available with raised side rails for additional patient safety.
All Mattresses come standard with a 4-way stretched nylon waterproof cover.

**48 HOUR
SHIPPING***

*order must be received by 2 PM EST

Mason Medical Products

Division of MRC Industries, Inc.

85 Denton Ave, New Hyde Park New York 11040
Tel: 800/233-4454 • 516/328-6900 • Fax: 516/328-6622




Multi-Ply™ Series 6500

VARIABLE DENSITY SLEEP SYSTEM

Multi-Ply 6500

Features and Benefits


- *Multi-layered* pressure reducing mattress with top layer of die-cut high density foam in head and torso area, and special high resilient foam in heel area, provides a comfortable pressure reducing environment
- Middle layer of foot section has precision channel cut foam which provides additional pressure reduction in heel area
- Bottom layer of high density foam provides a firm foundation
- Mason's SSPT Rail System™ provides solid-edge support, allowing easy patient transfer to and from bed
- Top Cover is a 4way stretch nylon which reduces heat buildup and provides a low shear non-friction surface that is anti-microbial and waterproof
- Border material & bottom panel of cover is 8.5 oz. grey Masongard™
- Designated head and foot sections
- Concealed zipper with Mason's *Barrier Stop™ Over-Flap* prevents liquids from contaminating mattress core
- Conforms to Cal. Tech. Bull. #117 for fire-retardancy



Multi-Ply 6500-Lite

Features and Benefits


- *Multi-layered* pressure reducing mattress with top layer of die-cut high density foam in the head and torso section for excellent pressure reduction
- Heal Section is multi-layered, the top layer consist of independent cell technology for superior pressure reduction
- Bottom layer of high density foam provides a firm foundation
- Mason's SSPT Rail System™ provides solid-edge support, allowing easy patient transfer to and from bed
- Top Cover is a 4way stretch nylon which reduces heat buildup and provides a low shear non-friction surface that is anti-microbial and waterproof
- Designated head and foot sections
- Concealed zipper with Mason's *Barrier Stop™ Over-Flap* prevents liquids from contaminating mattress core
- Conforms to Cal. Tech. Bull. #117 for fire-retardancy



Multi-Ply 6500-Global

Features and Benefits

- *Multi-layered* pressure reducing mattress that is soft enough for your lighter patients and durable enough for your heavier patients (450 lbs)
- The top layer of Lura-Quilt™ foam is aerated which allows air to continuously pass through reducing heat build up
- The sloped heel section reduces pressure and increases blood circulation by redistributing the load from the heel to the less vulnerable calf
- The heel section is precision die-cut to permit air to circulate under the patients heel, which keeps the skin cool and dry
- The Lura-Tex™ bottom layer is channel cut to allow continuous gatching of the head and foot section without complications
- Top Cover is a 4way stretch nylon which reduces heat buildup, low shear anti-microbial and water proof
- Designated head and foot sections
- Concealed zipper with Mason's *Barrier Stop™ Over-Flap* prevents liquids from contaminating mattress core
- Conforms to Cal. Tech. Bull. #117 for fire-retardancy



Mason Medical Products

Division of MRC Industries, Inc.

85 Denton Ave, New Hyde Park New York 11040
Tel: 800/233-4454 • 516/328-6900 • Fax: 516/328-6622

