

POWER SUPPORT SURFACE

EXCEL5000

LOW AIR LOSS/ALTERNATING PRESSURE OVERLAY SYSTEM



Mason Medical Products

Division of MRC Industries, Inc.

Tel: 800/233-4454 • 516/328-6900 • Fax: 516/328-6622

www.masonmedical.com

LOW AIR LOSS/ALTERNATING PRESSURE OVERLAY SYSTEM

The Excel 5000 combines low air loss and alternating pressure therapies for an extraordinary overlay system. It is an ideal system to prevent and treat pressure sores. The powerful, yet quiet air pump provides sophisticated, user-friendly operation. The Excel 5000 incorporates a system of 18 individual bladders. The pump offers an 8 minute cycle time; adjustable comfort settings for the individual patient's needs.



FEATURES AND BENEFITS

FOR THE HOME CARE AND LONG TERM MARKETS

MATTRESS & COVER

- 16 laser air holes to keep patient cool and dry.
- Removable, non-shear, anti-microbial fully enclosed cover encapsulates the entire 18 bladder mattress.
- CPR feature allows for rapid deflation of mattress system.
- Overall dimensions: 33" x 78.5" x 5"
- Weight: 9 lbs
- Alternating cell depth: 5"

EXCEL 5000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Soft/Firm Dial to obtain custom patient comfort pressure levels.
- Low pressure alarm (if mattress hose is disconnected).
- Adjustable hanging brackets.
- Hospital grade UL approved power cord.
- Kink resistant tubing.
- Weight Capacity : 300 lbs
- Voltage Range : 110-120/220-240 VAC
- Frequency Range : 50/60 Hz
- Dimension : 10" x 5" x 4"
- Weight : 3.5 lbs
- 12 Month Warranty

- # EX5000 EXCEL 5000 System
- # EX5120 EXCEL 5000 Pump (only)
- # EX5125 18 Bladder Mattress (only)
- # EX5140 Replacement Mattress Cover



Mason Medical Products
Division of MRC Industries, Inc.

Tel: 800/233-4454 • 516/328-6900 • Fax: 516/328-6622