



# Your Source for Quality Health Care Products



## MEDTRADE 2010 SHOW SPECIALS 5/12/10-5/13/10

### LS/AS 8800

LOW AIR LOSS & ALTERNATING PRESSURE MATTRESS



EO277

- Low air loss mattress with alternating pressure mode.
- Cell-on-cell design keeps bottom half of 8" bladders inflated in case of power failure.
- Static function to allow treatment of patient.
- Microprocessor pressure sensing technology to produce low pressure support.
- Pressure range adjusted by patients weight
- Warranty: 12 months.
- Weight capacity: 350lbs.

### AS 5000

ALTERNATING PRESSURE MATTRESS



EO277

- (20) 8" bladders for optimum pressure relief.
- Static function to allow treatment of patient.
- Soft/ firm dial to obtain custom comfort level by patients weight.
- Audio and visual low pressure alarms.
- Removable, non-shear, anti- microbial nylon cover.
- CPR valve for rapid deflation.
- Warranty: 12 months.
- Weight capacity: 350 lbs.

### LS 9500

TURN-STYLE LATERAL ROTATION WITH ON DEMAND LOW AIR MATTRESS



EO277

- On demand low air loss feature
- Lateral Rotation
- High airflow multi-function control unit
- 2" foam base in case of power failure
- 10, 20,30 & 60 minute adjustable turning cycles mode.
- 1/4, 1/2, 3/4 & full position adjustable turning positions
- Microprocessor pressure sensing technology to produce low pressure support.
- 10 adjustable comfort levels and adjustable turning angle up to 40°
- Audible and visual low pressure alarms.
- Vapor permeable, waterproof and fire retardant top cover.
- CPR deflation in 20 seconds.
- 12 month warranty.

### LS 9000B-42

BIATRIC LOW AIR LOSS/ALTERNATING PRESSURE MATTRESS



EO277

- Powerful Powerloft Pro Pump unit supports up to 1,000 lbs. (sizes 42", 48").
- Dual functions combining low air loss and alternating
- 6 airflow settings allow custom pressure selection for individual patient's needs.
- Static function to allow treatment of patient.
- Microprocessor pressure sensing technology produces low pressure
- 2 CPR valves for rapid deflation.
- Warranty: 12 months.

### SIERRA™ GEL OVERLAY



EO185

- Optimal environment for prevention/treatment of pressure sores
- One-piece construction requires no assembly
- Easy to clean, non shear/low friction surface
- Elastic corner straps to secure overlay to mattress
- Base Foundation is a high density solid one-piece foam constructed core
- Convuluted high density polyurethane foam topper

### 2000PP-II ALTERNATING PRESSURE OVERLAY SYSTEM



EO181

- Extremely durable, stand up unit
- Provides a maintenance free pressure reducing environment for the patient
- Lighted On/Off switch
- Air output: 4.2 liters
- Cycle Time: 5 minutes
- Power: 110V/60Hz

### BALANCED AIRE ADJUSTABLE AIR CUSHION



KO734

The Balanced Aire Adjustable air cushion is constructed of independent flexible interconnected air cells designed to immerse the patient to help disburse pressure evenly and increase blood flow to assist in the treatment and prevention of pressure sores. The waterproof, 4 way stretch cover provides low shear to help maintain sensitive skin tissue while the air pump will be used to adjust individual patient comfort. The non-skid bottom and adjustable strap offers additional security. Available in a low profile 2" cushion.

### PREMIER ONE™



E2601

Designed specifically to reduce interface pressure and maintain skin integrity. The entire cushion is made of Luratex™ foam. Luratex™ is a revolutionary patented high resilient aerated foam. The formulation allows air to pass through the foam to keep the patient cool, dry and comfortable. The covers are available in stretch knit or waterproof nylon

### PREMIER ONE™ BACK CUSHION



E2611

This top of line high resilient, high density contoured foam back cushion provides upper back comfort with advanced therapeutic lumbar support which will encourage correct patient posture and proper spine alignment. The cover is a 4 way stretch Nylon with a non skid back and comes complete with a solid insert to prevent hammocking along with 2 fully adjustable locking straps to secure cushion to the wheelchair

**Mason Medical Products**  
Division of MRC Industries, Inc.

85 Denton Ave, New Hyde Park, New York 11040  
Tel: 800/233-4454 • 516/328-6900 • Fax: 516/328-6622

