

A L T E R N A T I N G P R E S S U R E

Power Overlays

System 2000™
ALTERNATING PRESSURE SYSTEMS



2000PP-II
FIXED ALTERNATING
PRESSURE SYSTEM



2000PP-ADJ
ADJUSTABLE ALTERNATING
PRESSURE SYSTEM



Mason Medical Products
Division of MRC Industries, Inc.

Tel: 800/233-4454 • 516/328-6900 • Fax: 516/328-6622
www.masonmedical.com

ALTERNATING PRESSURE OVERLAYS

System 2000™ offers an economical approach to pressure sore management. The hospital grade alternating pressure pump and bubble pad is both user-friendly and maintenance free.



SYSTEM 2000PP-II FIXEDFLOW ALTERNATING PRESSURE PUMP

- Extremely durable, stand-up unit
- Provides a maintenance free pressure reducing environment for the patient
- Lighted On/Off switch • Air Output: 4.2 liters
- Cycle Time: 5 minutes • Power: 110V/60Hz.
- Power Cord: 10 foot w/ 3 pronged grounded plug
- UL® listed

BENEFITS AND FEATURES FOR THE HOME CARE AND LONG TERM MARKET

SYSTEM 2000™ BUBBLE PAD

Offers the patient more benefits than linear style pads.

The unique design of this bubble pad allows more touch points offering a superior therapeutic treatment.

- Material: Medical grade PVC with fire-retardant and bacteriostatic treatments
- Thickness: 11.5 gauge
- Design: Bubble Cells
- Size of Pad: 76"(L) x 31.5"(W) inflated
- Height of Cell: 2⁷/₈" inflated

- **2000PP-II 2002-II 2001** **System 2000 Fixed Flow System Fixed Flow Pump(only) Replacement Bubble Pads**

Universal Adaptor available for use with other brand name pumps



SYSTEM 2000PP-ADJ ADJUSTABLE AIR FLOW ALTERNATING PRESSURE PUMP

- Allows clinician to control air flow for patients weighing between 65 and 265 lbs.
- Manual adjustment dial for air loft of pad
- Hanging bed hooks • Lighted On/Off switch
- Air Output: 6 liters • Cycle Time: 5.5 minutes
- Power: 110V/60Hz.
- Power Cord: 10 foot with 3 pronged grounded plug
- UL® listed.

- **2000-ADJ 2002-ADJ 2001-FL** **System 2000 Adjustable System Adjustable Pump(only) Replacement Bubble Pads W/Flaps**